



Introducing the lalaPositioner



**SHONAQUIP
SOCIAL ENTERPRISE**
Enabling inclusion

Removable Foam Positioning Blocks



Durable Covers



Rigid Foldable Frame

lalaPositioner

The LalaPositioner is a useful positioning device designed to improve a child's postural alignment when spending time out of the wheelchair. It can be used to position the child while asleep during the night or for effective and comfortable positioning for short periods during the day. It helps to control severe muscle stiffness, spasms and uncontrolled movements and is particularly effective in preventing extreme flexed or extended positions of the body, head and limbs. It also allows for controlling the positioning of the head to improve the child's ability to breathe, swallow and socialise. Different users can be accommodated on the same positioner.

It is effective in promoting healthy growth and development of the hips, knees and spine in children by maintaining a normal and relaxed posture in lying, over longer periods.

LALAPOSTIONER SIZES & SPECS

	BABY	SMALL	MEDIUM	LARGE
Base & back cushions				
Base mattress size (width x length)	1000mm x 600mm	1200mm x 600mm	1500mm x 600mm	1800mm x 700mm
Base mattress thickness	75mm	75mm	100mm	100mm
Back mattress size (width x length)	1000mm x 350mm	1200mm x 350mm </td <td>1500mm x 450mm</td> <td>1800mm x 500mm</td>	1500mm x 450mm	1800mm x 500mm
Back mattress thickness	50mm	50mm	50mm	50mm
Posture support cushions				
Head cushion (L x W x H)	540 x 215 x 50mm	540 x 215 x 50mm	550 x 310 x 100mm	640 x 300 x 125mm
Leg Block Large (L x W x H)	400 x 300 x 50mm	550 x 350 x 75mm	600 x 385 x 100mm	700 x 390 x 150mm
Leg Block Small (L x W x H)	220 x 150 x 100mm	350 x 150 x 100mm	380 x 180 x 150mm	450 x 180 x 150mm

lala2 coming soon...



For any enquiries please contact our Cape Town Head Office

+27 21 797 8239 | info@ShonaquipSE.org.za | www.ShonaquipSE.org.za



**SHONAQUIP
SOCIAL ENTERPRISE**
Enabling inclusion