



jackStander

☐ Customised Device

Organisation: _____

Date: ____ DD ____ MM ____ YYYY

Client name: _____

Therapist: _____

Technician to call the therapist **BEFORE** making this device: No ☐ Yes ☐

Due Date: ____ DD ____ MM ____ YYYY

Key: ☐ Tickbox ☒ Upcharge Item


* The tables are a size reference guide for standard items.

BODY SUPPORT

	BABY ☐	SMALL ☐	MEDIUM ☐	LARGE ☐
Pelvic Width Range	160mm & 200mm	160mm & 200mm	210mm & 260mm	300mm & 350mm
Overall Length Range	750mm - 1100mm	1000mm - 1280mm	1100mm - 1450mm	1400mm - 1750mm
Heel to Pelvis Length	650mm	750mm	800mm	850mm
Tilt Angle Min & Max	90° - 180°	90° - 180°	90° - 180°	90° - 180°


OPTIONS

☒ **Optional:** Shoulder Side Support




STRAPS

	BABY	SMALL	MEDIUM	LARGE
Jack Stander Pelvic Strap	Size1	Size1	Size2	Size2
Chest Strap for Supine	Size1	Size2	Size2	Size3
Back Strap for Prone	Size1	Size1	Size1	Size2

Select: Jack Stander Pelvic Strap (W x L Strap)Size ☐ **1:** 100 x 1080mm ☐ **2:** 100 x 1310mm☐ **Optional:** Chest Strap for Supine positioning (Padding diameter)Size ☐ **1:** 100mm ☐ **2:** 120mm ☐ **3:** 140mm

Round chest strap

☐ **Optional:** Chest Strap Rectangular (W x L Padding)Size ☐ **1:** 120 x 260mm ☐ **2:** 130 x 300mm ☐ **3:** 130 x 300mm (Long Strap)☐ **Optional:** Back Strap for Prone positioning (W x H Padding)Size ☐ **1:** 200 x 150mm ☐ **2:** Custom

Prone back strap

☐ **Custom:** strap

**Remember to indicate:

- Webbing widths and lengths.
- Foam padding type, widths, lengths and thickness.
- Types of buckles to use.
- Draw a diagram where possible.

SPECIAL INSTRUCTIONS

